

Boule Con-5Diff A1

CONTROL

IVD

Tri-Level	72605-1K	2026/09/05
Low	72605-11	2026/09/05
Normal	72605-12	2026/09/05
High	72605-13	2026/09/05

Open vial stability: 14 days ¹

Medonic M51 / Swelab Lumi

Parameter ²	Unit ³	Low		Normal		High	
		Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵
WBC	10 ⁹ /L	3,50	± 0.50	8,50	± 1.00	19,97	± 2.50
Neu%	%	48,9	16,0	59,3	8,0	66,2	8,0
Lym%	%	40,4	21,0	28,2	8,0	20,7	12,0
Mon%	%	5,0	5,0	5,3	5,0	4,7	6,0
Eos%	%	5,7	5,0	7,2	6,0	8,4	8,0
Bas%	%	0,7	10,0	2,8	10,0	2,4	10,0
Neu#	10 ⁹ /L	1,73	0,50	5,14	0,70	13,46	2,00
Lym#	10 ⁹ /L	1,43	0,70	2,45	0,80	4,22	2,00
Mon#	10 ⁹ /L	0,18	0,30	0,46	0,50	0,95	1,10
Eos#	10 ⁹ /L	0,20	0,20	0,62	0,60	1,72	1,40
Bas#	10 ⁹ /L	0,03	0,40	0,24	0,80	0,48	1,80
RBC	10 ¹² /L	2,30	0,24	4,61	0,30	5,30	0,50
HGB	g/L	65	6	139	6	177	8
HCT	%	19,8	2,0	42,2	3,0	52,3	4,0
MCV	fL	86,2	5,0	91,6	5,0	98,7	6,0
MCH	pg	28,3	2,5	30,2	2,5	33,4	2,5
MCHC	g/L	328	30	329	30	338	30
RDW-CV	%	16,2	3,0	15,0	3,0	14,5	3,0
PLT	10 ⁹ /L	61	20	268	40	486	60
MPV	fL	8,0	3,0	7,3	3,0	7,8	3,0

No mean values and ranges are provided for research parameters-RDW-SD, PDW, PCT, P-LCR, and P-LCC. User may establish these based on their usage.



<p>1 bg Трайност на отворен флакон: 14 дни cs Stabilita otevřené ampulky: 14 dní da Åben rør stabilitet: 14 dage de Stabilität geöffneter Flaschen: 14 Tage el Διατήρηση δείγματος μετά την αποσφράγιση: 14 ημέρα es Estabilidad de la cápsula abierta: 14 días et Stabiilsus avatud vialis: 14 päevade arv fr Stabilité en flacon ouvert: 14 jour hr Stabilnost otvorene bočice: 14 dani it Stabilità della fiala aperta: 14 giorni lt Stabilumas atidarius buteliukoje: 14 dienų lv Atvērt flakonū stabilitāte: 14 dienas no Åpen rør/glass stabilitet: 14 dager pl Trwałość po otwarciu fiołki: 14 Liczba dni pt Estabilidade após abertura do frasco: 14 dias ro Termenul de valabilitate al fiolei desfăcute: 14 zile ru Стабильность вскрытого флакона: 14 дней sk Stabilita otvorenej injekčnej liekovky: 14 dní sl Stabilnost odprte stekleničke: 14 dni sr Proizvod stabilan nakon otvaranja pakovanja: 14 dana sv Hållbarhet för öppnad flaska: 14 dagar tr Açık şişe stabilitesi gün: 14 gün</p>	<p>2 Параметър Parameter Parameter Parameter Παράμετρος Parámetro Parameeter Paramètre Parametar Parametro Parametras Parameters Parameter Kryterium Parâmetro Parametru Параметр Parameter Parameter Parameter Parameter Parameter Parameter</p>	<p>3 Единица Jednotka Enhed Einheit Μονάδα Unidad Ühik Unité Jedinica Unità Vienetas Vienība Enhet Jednostka Unidade Unitate Единица измерения Jednotka Enota Jedinica Enhet Birim</p>	<p>4 Средно Průměr Middelværdi Mittelwert Μέση Media Keskmine Moyenne Sredstvo Media Vidurkis Nozīmē Middelverdi Średnia Média Medie Среднее Priemer Povprečno Srednja vrednost Medelvärde Ortalama</p>	<p>5 Обхват Rozsah Område Bereich Μέση Αναμενόμενο εύρος α Intervalo Vahemik Intervalle Raspon Intervallo Intervalas Diapazons Område Zakres Intervalo Interval диапазон Rozsah Območje Opseg Intervall Aralığı</p>
--	--	---	--	--



www.boule.com

Boule
Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

31340-14