

LOT: 12605-14
EXP: 2026-08-04
Open Vial Stability 5 days

EXIGO

| | WBC 10 ⁹ /L | RBC 10 ¹² /L | Hgb g/dL | MCV fL | PLT 10 ⁹ /L |
|------------------------|------------------------|-------------------------|----------|--------|------------------------|
| System specific values | 8.5 | 4.21 | 11.0 | 78.8 | 223 |
| Tolerance limits ± | 0.3 | 0.10 | 0.2 | 2.0 | 15 |



For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12605-14
Prod.: Vet-Cal
ID: 1260514+
Exp.: 4-Aug-2026

