

LOT: 12605-52
EXP: 2026-10-18
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	10 ⁹ /L	8.6	0.6	8.0	9.2
Hct	%	33.2	2.5	30.7	35.7
MCV	fL	79.2	5.0	74.2	84.2
RBC	10 ¹² /L	4.19	0.18	4.01	4.37
Hgb	g/dL	11.0	0.4	10.6	11.4
MCH	pg	26.3	1.5	24.8	27.8
MCHC	g/dL	33.1	2.2	30.9	35.3
PLT	10 ⁹ /L	221	30	191	251



For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12605-52
Prod.: Vet-Con Norm
ID: 1260552+
Exp.: 18-Oct-2026

