

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22606-1K	2026-11-02
Low	22606-11	2026-11-02
Normal	22606-12	2026-11-02
High	22606-13	2026-11-04

MEDONIC

Medonic M-series/M32

CONTROL		L	
LOT		22606-11	
		Mean ²	± Range ³
RBC	10 ¹² /L	2.22	0.12
MCV	fL	69.2	5.0
RDW	%	14.9	4.0
Hct	%	15.4	2.5
Plt	10 ⁹ /L	82	15
MPV	fL	11.2	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	25.2	1.5
MCHC	g/dL	36.4	3.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	43.3	7.0
Mid	%	5.7	5.0
Gran	%	51.0	7.0

CONTROL		N	
LOT		22606-12	
		Mean ²	± Range ³
		4.10	0.18
		84.1	5.0
		13.0	4.0
		34.5	3.0
		240	30
		11.2	1.5
		8.5	0.6
		12.0	0.4
		29.3	1.5
		34.8	3.0
		4.0	0.6
		0.5	0.5
		4.0	1.0
		46.7	5.0
		6.2	5.0
		47.1	5.0

CONTROL		H	
LOT		22606-13	
		Mean ²	± Range ³
		5.00	0.22
		93.5	5.0
		13.2	4.0
		46.8	3.5
		528	60
		11.3	1.5
		21.0	1.8
		15.9	0.5
		31.8	1.5
		34.0	3.0
		12.8	1.5
		1.3	1.0
		6.9	2.5
		61.3	5.0
		6.0	5.0
		32.7	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dni (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fioleki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelverdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı



www.boule.com

Boule



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25121-8

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series
Lot: 22606-11
Prod.: Con-Diff Low
ID: 2260611+
Exp.: 2-Nov-2026

Instr.: M-series
Lot: 22606-12
Prod.: Con-Diff Norm
ID: 2260612+
Exp.: 2-Nov-2026

Instr.: M-series
Lot: 22606-13
Prod.: Con-Diff High
ID: 2260613+
Exp.: 4-Nov-2026



Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22606-1K	2026-11-02
Low	22606-11	2026-11-02
Normal	22606-12	2026-11-02
High	22606-13	2026-11-04

SWELAB

		CONTROL	L
		LOT	22606-11
		Mean ²	± Range ³
RBC	10 ¹² /L	2.26	0.12
MCV	fL	70.1	5.0
RDW	%	14.0	4.0
Hct	%	15.8	2.5
Plt	10 ⁹ /L	83	15
MPV	fL	11.8	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.7	0.3
MCH	pg	25.2	1.5
MCHC	g/dL	36.1	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	47.3	7.0
Mid	%	3.7	5.0
Gran	%	49.0	7.0

		CONTROL	N
		LOT	22606-12
		Mean ²	± Range ³
		4.11	0.18
		85.2	5.0
		11.7	4.0
		35.0	3.0
		242	30
		11.8	1.5
		8.6	0.6
		12.0	0.4
		29.2	1.5
		34.3	3.0
		4.3	0.6
		0.4	0.5
		3.9	1.0
		50.3	5.0
		4.2	5.0
		45.5	5.0

		CONTROL	H
		LOT	22606-13
		Mean ²	± Range ³
		5.02	0.22
		94.8	5.0
		12.1	4.0
		47.6	3.5
		551	60
		11.9	1.5
		21.1	1.8
		16.0	0.5
		31.9	1.5
		33.6	3.0
		13.7	1.5
		0.8	1.0
		6.6	2.5
		64.6	5.0
		3.9	5.0
		31.5	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dni (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilitumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiolek 14 dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcută 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Sredno
 - Sredstvo
 - Průměr
 - Middelvaerd
 - Keskmine
 - Moyenne
 - Mittelwert
 - Mέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı



www.boule.com



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25121-8

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
Lot: 22606-11
Prod.: Con-Diff Low
ID: 2260611+
Exp.: 2-Nov-2026

Instr.: Swelab Alfa
Lot: 22606-12
Prod.: Con-Diff Norm
ID: 2260612+
Exp.: 2-Nov-2026

Instr.: Swelab Alfa
Lot: 22606-13
Prod.: Con-Diff High
ID: 2260613+
Exp.: 4-Nov-2026

