

IVD

5 days ¹

Quintus							
Parameter ²	RBC 10 ¹² /L	MCV fL	Plt 10 ⁹ /L	MPV fL	WBC 10 ⁹ /L	Hgb g/dL	RDW _{cv} %
Mean ³	4.08	84.8	233	9.5	8.6	13.1	14.5
Range ± ⁴	0.10	2.0	20	1.5	0.3	0.3	1.5

Target X	Target Y
87	116

Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Αράλδι



www.boule.com



Boule



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

26577-6

For Quintus with barcode reader

1. In the Main Menu go to “*Calibration QC*” / “*Enter QC reference*”.
2. Highlight the first Barcode field, scan in each calibrator/control barcode of the column in sequence (1-9) and click “*OK*”.
3. Repeat steps 1 and 2 for each lot to be entered.

Instr.: Quintus
Lot: 62607-44
Prod.: Calibrator
ID: 6260744+
Exp.: 2-Oct-2026



4.