

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22608-1K	2027-01-04
Low	22608-11	2027-01-04
Normal	22608-12	2027-01-04
High	22608-13	2027-01-04

MEDONIC

CONTROL		L	
LOT		22608-11	
Medonic M-series/M32		Mean ²	± Range ³
RBC	10 ¹² /L	2.25	0.12
MCV	fL	70.7	5.0
RDW	%	14.1	4.0
Hct	%	15.9	2.5
Plt	10 ⁹ /L	83	15
MPV	fL	11.3	1.5
WBC	10 ⁹ /L	3.2	0.3
Hgb	g/dL	5.8	0.3
MCH	pg	25.8	1.5
MCHC	g/dL	36.5	3.0
Lymph#	10 ⁹ /L	1.4	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.6	0.8
Lymph	%	43.5	7.0
Mid	%	5.7	5.0
Gran	%	50.8	7.0

CONTROL		N	
LOT		22608-12	
		Mean ²	± Range ³
		4.14	0.18
		82.6	5.0
		13.4	4.0
		34.2	3.0
		231	30
		11.7	1.5
		8.6	0.6
		12.1	0.4
		29.2	1.5
		35.4	3.0
		4.0	0.6
		0.5	0.5
		4.1	1.0
		46.2	5.0
		6.3	5.0
		47.5	5.0

CONTROL		H	
LOT		22608-13	
		Mean ²	± Range ³
		4.98	0.22
		95.3	5.0
		12.9	4.0
		47.5	3.5
		548	60
		12.0	1.5
		21.5	1.8
		16.3	0.5
		32.7	1.5
		34.3	3.0
		13.2	1.5
		1.2	1.0
		7.1	2.5
		61.3	5.0
		5.8	5.0
		32.9	5.0

1

bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilnost otvorene bočice 14 dani (hr)
cs	Stabilita otevřené ampulky 14 dny (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud vialis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διάρθρωση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiolek 14 Liczba dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvoreneje injekčneje liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 dias (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2

Sredno
Sredstvo
Průměr
Middelværdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nozīmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3

Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Aralığı



www.boule.com



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25121-8

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series

Lot: 22608-11

Prod.: Con-Diff Low

ID: 2260811+

Exp.: 4-Jan-2027

Instr.: M-series

Lot: 22608-12

Prod.: Con-Diff Norm

ID: 2260812+

Exp.: 4-Jan-2027

Instr.: M-series

Lot: 22608-13

Prod.: Con-Diff High

ID: 2260813+

Exp.: 4-Jan-2027

1:



Boule CBA V2



Boule CBA V2



Boule CBA V2

2:



2260811+



2260812+



2260813+

3:



27010400001153030000



27010400002153030000



27010400003153030000

4:



10004029300000000000



10004029800000000000



10004029700000000000

5:



22512707500580315925



41418826501210434230



49822953501630547535

6:



08315113152581536530



23130117152921535430



54860120153271534330

7:



03203014040020301608



08606040060050504110



21518132150121007125

8:



14140435700575050870



13440462500635047550



12940613500585032950

9:



11116926914938375892



11111629083474513056



11111759936886062773

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22608-1K	2027-01-04
Low	22608-11	2027-01-04
Normal	22608-12	2027-01-04
High	22608-13	2027-01-04

SWELAB

CONTROL		L	
LOT		22608-11	
Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.26	0.12
MCV	fL	71.8	5.0
RDW	%	13.2	4.0
Hct	%	16.2	2.5
Plt	10 ⁹ /L	82	15
MPV	fL	11.8	1.5
WBC	10 ⁹ /L	3.2	0.3
Hgb	g/dL	5.9	0.3
MCH	pg	26.1	1.5
MCHC	g/dL	36.4	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.5	0.8
Lymph	%	48.0	7.0
Mid	%	3.8	5.0
Gran	%	48.2	7.0

CONTROL		N	
LOT		22608-12	
		Mean ²	± Range ³
		4.15	0.18
		84.0	5.0
		12.2	4.0
		34.9	3.0
		234	30
		12.4	1.5
		8.6	0.6
		12.1	0.4
		29.2	1.5
		34.7	3.0
		4.3	0.6
		0.4	0.5
		3.9	1.0
		50.4	5.0
		4.1	5.0
		45.5	5.0

CONTROL		H	
LOT		22608-13	
		Mean ²	± Range ³
		4.93	0.22
		96.7	5.0
		11.9	4.0
		47.7	3.5
		558	60
		12.7	1.5
		21.2	1.8
		16.3	0.5
		33.1	1.5
		34.2	3.0
		13.7	1.5
		0.8	1.0
		6.7	2.5
		64.5	5.0
		3.7	5.0
		31.8	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
- hr Stabilitet otvorene bočice 14 dni (hr)
- cs Stabilita otevřených ampulky 14 dní (cs)
- da Åben rør stabilitet 14 dage (da)
- et Stabiilsus avatud vialis 14 päevade (et)
- fr Stabilité en flacon ouvert 14 jour (fr)
- de Stabilität geöffneter Flaschen 14 tage (de)
- el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
- it Stabilità della fiala aperta 14 giorni (it)
- lv Atvērt flakonu stabilitāti 14 dienas (lv)
- lt Stabilitumas atidarius buteliuką 14 dienas (lt)
- no Åpen rør/glass stabilitet 14 dager (no)
- pl Trwałość po otwarciu fiolek 14 dni (pl)
- pt Estabilidade após abertura do frasco 14 dias (pt)
- ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
- ru Стабильность вскрытого флакона 14 дней (ru)
- sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
- sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
- sl Stabilita otvoreneje injekčneje liekovky 14 dni (sl)
- es Estabilidad de la cápsula abierta 14 días (es)
- sv Hållbarhet för öppnad flaska 14 dagar (sv)
- tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
- Sredstvo
- Průměr
- Middelverdi
- Keskmine
- Moyenne
- Mittelwert
- Μέση
- Media
- Nožimě
- Vidurkis
- Middelverdi
- Średnia
- Média
- Medie
- Среднее
- Srednja vrednost
- Priemer
- Povprečno
- Media
- Medelvärde
- Ortalama

- 3
- Обхват
- Raspon
- Rozsah
- Område
- Vahemik
- Intervalle
- Bereich
- Αναμενόμενο εύρος α
- Intervallo
- Diapazons
- Intervalas
- Område
- Zakres
- Intervalo
- Interval
- диапазон
- Opseg
- Rozsah
- Območje
- Intervalo
- Intervall
- Aralığı



www.boule.com



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25121-8

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa

Lot: 22608-11

Prod.: Con-Diff Low

ID: 2260811+

Exp.: 4-Jan-2027

Instr.: Swelab Alfa

Lot: 22608-12

Prod.: Con-Diff Norm

ID: 2260812+

Exp.: 4-Jan-2027

Instr.: Swelab Alfa

Lot: 22608-13

Prod.: Con-Diff High

ID: 2260813+

Exp.: 4-Jan-2027

