

LOT: 12608-52
EXP: 2027-01-17
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	$10^9/L$	8.5	0.6	7.9	9.1
Hct	%	31.7	2.5	29.2	34.2
MCV	fL	77.1	5.0	72.1	82.1
RBC	$10^{12}/L$	4.11	0.18	3.93	4.29
Hgb	g/dL	10.6	0.4	10.2	11.0
MCH	pg	25.8	1.5	24.3	27.3
MCHC	g/dL	33.4	2.2	31.2	35.6
PLT	$10^9/L$	223	30	193	253



For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12608-52
Prod.: Vet-Con Norm
ID: 1260852+
Exp.: 17-Jan-2027

