

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22609-1K	2027-02-01
Low	22609-11	2027-02-01
Normal	22609-12	2027-02-01
High	22609-13	2027-02-03

MEDONIC

CONTROL		L	
LOT		22609-11	
Medonic M-series/M32		Mean ²	± Range ³
RBC	10 ¹² /L	2.23	0.12
MCV	fL	66.6	5.0
RDW	%	14.9	4.0
Hct	%	14.9	2.5
Plt	10 ⁹ /L	84	15
MPV	fL	11.7	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.4	0.3
MCH	pg	24.2	1.5
MCHC	g/dL	36.2	3.0
Lymph#	10 ⁹ /L	1.4	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	43.7	7.0
Mid	%	5.4	5.0
Gran	%	50.9	7.0

CONTROL		N	
LOT		22609-12	
		Mean ²	± Range ³
		4.17	0.18
		82.6	5.0
		13.3	4.0
		34.4	3.0
		234	30
		11.7	1.5
		8.7	0.6
		12.1	0.4
		29.0	1.5
		35.2	3.0
		4.0	0.6
		0.6	0.5
		4.1	1.0
		46.7	5.0
		6.5	5.0
		46.8	5.0

CONTROL		H	
LOT		22609-13	
		Mean ²	± Range ³
		5.07	0.22
		92.8	5.0
		13.0	4.0
		47.0	3.5
		513	60
		11.8	1.5
		20.8	1.8
		16.2	0.5
		32.0	1.5
		34.5	3.0
		12.8	1.5
		1.2	1.0
		6.8	2.5
		61.4	5.0
		5.8	5.0
		32.8	5.0

1

bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilnost otvorene bočice 14 dni (hr)
cs	Stabilita otevřené ampulky 14 dny (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud vialis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διάρθρωση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiolek 14 Liczba dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvorene injekčneje liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 dias (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2

Средно
Sredstvo
Průměr
Middelværdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nozīmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3

Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Aralığı



www.boule.com

Boule



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25121-8

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series

Lot: 22609-11

Prod.: Con-Diff Low

ID: 2260911+

Exp.: 1-Feb-2027

Instr.: M-series

Lot: 22609-12

Prod.: Con-Diff Norm

ID: 2260912+

Exp.: 1-Feb-2027

Instr.: M-series

Lot: 22609-13

Prod.: Con-Diff High

ID: 2260913+

Exp.: 3-Feb-2027

1:



Boule CBA V2



Boule CBA V2



Boule CBA V2

2:



2260911+



2260912+



2260913+

3:



27020100001153030000



27020100002153030000



27020300003153030000

4:



10004029500000000000



10004029800000000000



10004029700000000000

5:



22312666500540314925



41718826501210434430



50722928501620547035

6:



08415117152421536230



23430117152901535230



51360118153201534530

7:



03303014040020301708



08706040060060504110



20818128150121006825

8:



14940437700545050970



13340467500655046850



13040614500585032850

9:



11117881508608947234



11111376447393671812



11118224366748974148

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22609-1K	2027-02-01
Low	22609-11	2027-02-01
Normal	22609-12	2027-02-01
High	22609-13	2027-02-03

SWELAB

CONTROL		L	
LOT		22609-11	
Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.23	0.12
MCV	fL	67.7	5.0
RDW	%	14.0	4.0
Hct	%	15.1	2.5
Plt	10 ⁹ /L	85	15
MPV	fL	12.3	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	24.7	1.5
MCHC	g/dL	36.4	3.0
Lymph#	10 ⁹ /L	1.7	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.6	0.8
Lymph	%	47.6	7.0
Mid	%	3.9	5.0
Gran	%	48.5	7.0

CONTROL		N	
LOT		22609-12	
		Mean ²	± Range ³
		4.13	0.18
		83.4	5.0
		12.1	4.0
		34.4	3.0
		242	30
		12.4	1.5
		8.5	0.6
		12.0	0.4
		29.1	1.5
		34.9	3.0
		4.3	0.6
		0.3	0.5
		3.9	1.0
		50.6	5.0
		4.1	5.0
		45.3	5.0

CONTROL		H	
LOT		22609-13	
		Mean ²	± Range ³
		5.04	0.22
		94.1	5.0
		12.0	4.0
		47.4	3.5
		524	60
		12.5	1.5
		20.6	1.8
		16.2	0.5
		32.1	1.5
		34.2	3.0
		13.3	1.5
		0.7	1.0
		6.6	2.5
		64.5	5.0
		3.6	5.0
		31.9	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
- hr Stabilität otvorene bočice 14 dni (hr)
- cs Stabilita otevřených ampulky 14 dní (cs)
- da Åben rør stabilitet 14 dage (da)
- et Stabiilsus avatud viaalis 14 päevade (et)
- fr Stabilité en flacon ouvert 14 jour (fr)
- de Stabilität geöffneter Flaschen 14 tage (de)
- el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
- it Stabilità della fiala aperta 14 giorni (it)
- lv Atvērt flakonu stabilitāti 14 dienas (lv)
- lt Stabilitumas atidarius buteliuką 14 dienas (lt)
- no Åpen rør/glass stabilitet 14 dager (no)
- pl Trwałość po otwarciu fiolek 14 dni (pl)
- pt Estabilidade após abertura do frasco 14 dias (pt)
- ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
- ru Стабильность вскрытого флакона 14 дней (ru)
- sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
- sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
- sl Stabilita otvoreneje injekčneje liekovky 14 dni (sl)
- es Estabilidad de la cápsula abierta 14 días (es)
- sv Hållbarhet för öppnad flaska 14 dagar (sv)
- tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
- Sredstvo
- Průměr
- Middelverdi
- Keskmine
- Moyenne
- Mittelwert
- Μέση
- Media
- Nožimě
- Vidurkis
- Middelverdi
- Średnia
- Média
- Medie
- Среднее
- Srednja vrednost
- Priemer
- Povprečno
- Media
- Medelvärde
- Ortalama

- 3
- Обхват
- Raspon
- Rozsah
- Område
- Vahemik
- Intervalle
- Bereich
- Αναμενόμενο εύρος α
- Intervallo
- Diapazons
- Intervalas
- Område
- Zakres
- Intervalo
- Interval
- диапазон
- Opseg
- Rozsah
- Območje
- Intervalo
- Intervall
- Aralığı



www.boule.com

Boule



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25121-8

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa

Lot: 22609-11

Prod.: Con-Diff Low

ID: 2260911+

Exp.: 1-Feb-2027

Instr.: Swelab Alfa

Lot: 22609-12

Prod.: Con-Diff Norm

ID: 2260912+

Exp.: 1-Feb-2027

Instr.: Swelab Alfa

Lot: 22609-13

Prod.: Con-Diff High

ID: 2260913+

Exp.: 3-Feb-2027

1:



B o u l e C B A V 3



B o u l e C B A V 3



B o u l e C B A V 3

2:



2260911+



2260912+



2260913+

3:



27020100001153030000



27020100002153030000



27020300003153030000

4:



20014018000000000000



20014018000000000000



20014018000000000000

5:



22312677500550315125



41318834501200434430



50422941501620547435

6:



08515123152471536430



24230124152911534930



52460125153211534230

7:



03403017040010301608



08506043060030503910



20618133150071006625

8:



14040476700395048570



12140506500415045350



12040645500365031950

9:



11117621263914131311



11116699318130619413



11110286235658489597